

Calcium Sensing Receptor protein (STJ503814)

STJ503814

GENERAL INFORMATION

Product Type	Proteins
Short Description	Calcium Sensing Receptor-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	846
Gene Symbol	CASR
Uniprot ID	CASR_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	