

CD36 protein (STJ503833)
STJ503833**GENERAL INFORMATION**

Product Type	Proteins
Short Description	CD36-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	948
Gene Symbol	CD36
Uniprot ID	CD36_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human CD36 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	