

CD4 protein (STJ503834)

STJ503834

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	920
Gene Symbol	CD4
Uniprot ID	CD4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 270-320 on human CD4 protein.
Immunogen Region	270-320
Specificity	
Immunogen Sequence	