

CD86 protein (STJ503837)

STJ503837

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD86-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	942
Gene Symbol	CD86
Uniprot ID	CD86_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-329 on human CD86 protein.
Immunogen Region	250-329
Specificity	
Immunogen Sequence	