

CD99 protein (STJ503838)
STJ503838**GENERAL INFORMATION**

Product Type	Proteins
Short Description	CD99-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 50-100 to human CD99 antigen.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	