

CDX2 protein (STJ503841)

STJ503841

GENERAL INFORMATION

Product Type	Proteins
Short Description	CDX2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1045
Gene Symbol	CDX2
Uniprot ID	CDX2_HUMAN
Immunogen	Mid-Region epitope synthetic peptide taken within amino acid region 100-200 on human CDX2 protein.
Immunogen Region	100-200
Specificity	
Immunogen Sequence	