

Cofilin protein (STJ503842)

STJ503842

GENERAL INFORMATION

Product Type	Proteins
Short Description	Cofilin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1072
Gene Symbol	CFL1
Uniprot ID	COF1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 120-170 on human Cofilin protein.
Immunogen Region	120-170
Specificity	
Immunogen Sequence	