

C6orf150 protein (STJ503843)

STJ503843

GENERAL INFORMATION

Product Type	Proteins
Short Description	C6orf150-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	214763
Gene Symbol	Cgas
Uniprot ID	CGAS_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 450-500 on mouse cGAS protein.
Immunogen Region	450-500
Specificity	
Immunogen Sequence	