

Chimaerin beta 3 protein (STJ503847)

STJ503847

GENERAL INFORMATION

Product Type	Proteins
Short Description	Chimaerin beta 3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1124
Gene Symbol	CHN2
Uniprot ID	CHIO_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-100 on human beta-chimaerin isoform 3 protein.
Immunogen Region	1-100
Specificity	
Immunogen Sequence	