

CHERP protein (STJ503848)

STJ503848

GENERAL INFORMATION

Product Type	Proteins
Short Description	CHERP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10523
Gene Symbol	CHERP
Uniprot ID	CHERP_HUMAN
Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 846-913 of human CHERP. Sequence: EQGIQDPIKG GDVRDKWDQY KGVGVALLDDP YENYRRNKSY SFIARMKARD
Immunogen Region	846-913
Specificity	
Immunogen Sequence	