

Cytokeratin 20 protein (STJ503850)

STJ503850

GENERAL INFORMATION

Product Type	Proteins
Short Description	Cytokeratin 20-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54474
Gene Symbol	KRT20
Uniprot ID	K1C20_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 human keratin, type I cytoskeletal 20 protein.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	