

cKit protein (STJ503859)
STJ503859**GENERAL INFORMATION**

Product Type	Proteins
Short Description	cKit-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3815
Gene Symbol	KIT
Uniprot ID	KIT_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 915-965 on human Tyrosine-Kinase c-Kit protein.
Immunogen Region	915-965
Specificity	
Immunogen Sequence	