

Calcitonin protein (STJ503863)

STJ503863

GENERAL INFORMATION

Product Type	Proteins
Short Description	Calcitonin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	796
Gene Symbol	CALCA
Uniprot ID	CALC_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 30-80 on human Calcitonin protein.
Immunogen Region	30-80
Specificity	
Immunogen Sequence	