

Cornichon protein (STJ503869)

STJ503869

GENERAL INFORMATION

Product Type	Proteins
Short Description	Cornichon-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10175
Gene Symbol	CNIH1
Uniprot ID	CNIH1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 90-120 on human Cornichon protein.
Immunogen Region	90-120
Specificity	
Immunogen Sequence	