

GPIP137 protein (STJ503870)

STJ503870

GENERAL INFORMATION

Product Type	Proteins
Short Description	GPIP137-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4076
Gene Symbol	CAPRIN1
Uniprot ID	CAPR1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-250 on human Caprin 1 protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	