

CTLA4 protein (STJ503874)

STJ503874

GENERAL INFORMATION

Product Type	Proteins
Short Description	CTLA4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1493
Gene Symbol	CTLA4
Uniprot ID	CTLA4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 150-200 on human CTLA4 protein
Immunogen Region	150-200
Specificity	
Immunogen Sequence	