

beta Catenin protein (STJ503875)
STJ503875**GENERAL INFORMATION**

Product Type	Proteins
Short Description	beta Catenin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1499
Gene Symbol	CTNNB1
Uniprot ID	CTNB1_HUMAN
Immunogen	N-epitope synthetic peptide taken within amino acid region 50-100 on human beta CTNNB protein.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	