

CTR1/SLC31A1 protein (STJ503877)

STJ503877

GENERAL INFORMATION

Product Type	Proteins
Short Description	CTR1/SLC31A1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1317
Gene Symbol	SLC31A1
Uniprot ID	COPT1_HUMAN
Immunogen	N-epitope synthetic peptide taken within amino acid region 1-100 on human High affinity copper uptake protein 1.
Immunogen Region	1-100
Specificity	
Immunogen Sequence	