

CXCR1 protein (STJ503880)

STJ503880

GENERAL INFORMATION

Product Type	Proteins
Short Description	CXCR1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	227288
Gene Symbol	Cxcr1
Uniprot ID	CXCR1_MOUSE
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on mouse CXCR1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	