

CXCR2 protein (STJ503881)

STJ503881

GENERAL INFORMATION

Product Type	Proteins
Short Description	CXCR2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3579
Gene Symbol	CXCR2
Uniprot ID	CXCR2_HUMAN
Immunogen	Synthetic peptide: EDLSNYSYSTLPPFLDAAAPC Corresponding to amino acids 18-39 of Human CXCR2.
Immunogen Region	18-39
Specificity	
Immunogen Sequence	