

CXCR2 protein (STJ503882)

STJ503882

GENERAL INFORMATION

Product Type	Proteins
Short Description	CXCR2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	729230
Gene Symbol	CCR2
Uniprot ID	CCR2_HUMAN
Immunogen	C-epitope Synthetic peptide corresponding to unique amino acid sequences on IL8 Receptor beta protein.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	