

CXCR4 protein (STJ503884)

STJ503884

GENERAL INFORMATION

Product Type	Proteins
Short Description	CXCR4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7852
Gene Symbol	CXCR4
Uniprot ID	CXCR4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 15-65 on CXCR4 protein.
Immunogen Region	15-65
Specificity	
Immunogen Sequence	