

Cyclin D1 protein (STJ503892)

STJ503892

GENERAL INFORMATION

Product Type	Proteins
Short Description	Cyclin D1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	595
Gene Symbol	CCND1
Uniprot ID	CCND1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acid sequences 268-295 on human Cyclin D1 protein.
Immunogen Region	268-295
Specificity	
Immunogen Sequence	