

CTR1/SLC31A1 protein (STJ503925)

STJ503925

GENERAL INFORMATION

Product Type	Proteins
Short Description	CTR1/SLC31A1-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1317
Gene Symbol	SLC31A1
Uniprot ID	COPT1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	