

CXCR7 protein (STJ503932)
STJ503932**GENERAL INFORMATION**

Product Type	Proteins
Short Description	CXCR7-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57007
Gene Symbol	ACKR3
Uniprot ID	ACKR3_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	