

C6orf150 protein (STJ503935)

STJ503935

GENERAL INFORMATION

Product Type	Proteins
Short Description	C6orf150-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	115004
Gene Symbol	CGAS
Uniprot ID	CGAS_HUMAN
Immunogen	Synthetic peptide taken within amino acid regions 350-450 on human cGAS protein.
Immunogen Region	350-450
Specificity	
Immunogen Sequence	