

## Dopamine Receptor D1 protein (STJ503938)

STJ503938

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Dopamine Receptor D1-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	NA
<b>Gene Symbol</b>	Drd1
<b>Uniprot ID</b>	DRD1_RAT
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 250-300 on mouse D (1A) dopamine receptor.
<b>Immunogen Region</b>	250-300
<b>Specificity</b>	
<b>Immunogen Sequence</b>	