

Dopamine Receptor D1 protein (STJ503940)

STJ503940

GENERAL INFORMATION

Product Type	Proteins
Short Description	Dopamine Receptor D1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Drd1
Uniprot ID	DRD1_RAT
Immunogen	Synthetic peptide taken within amino acid region 290-340 on human D-1 dopamine receptor protein.
Immunogen Region	290-340
Specificity	
Immunogen Sequence	