

Dopamine D2 Receptor protein (STJ503941)

STJ503941

GENERAL INFORMATION

Product Type	Proteins
Short Description	Dopamine D2 Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 150-200 on rat Dopamine Receptor D2 protein.
Immunogen Region	150-200
Specificity	
Immunogen Sequence	