

Dopamine D2 Receptor protein (STJ503942)

STJ503942

GENERAL INFORMATION

Product Type	Proteins
Short Description	Dopamine D2 Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1813
Gene Symbol	DRD2
Uniprot ID	DRD2_HUMAN
Immunogen	Synthetic peptide corresponding to unique sequence on human Dopamine D2 Receptor protein.
Immunogen Region	
Specificity	
Immunogen Sequence	