

Dopamine Receptor D3 protein (STJ503943)

STJ503943

GENERAL INFORMATION

Product Type	Proteins
Short Description	Dopamine Receptor D3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	29238
Gene Symbol	Drd3
Uniprot ID	DRD3_RAT
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on Dopamine Receptor D3 protein, peptide was covalently modified. Sequence: CQACHVSPELYRATTWGY
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	