

Dopamine D5 Receptor protein (STJ503945)

STJ503945

GENERAL INFORMATION

Product Type	Proteins
Short Description	Dopamine D5 Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13492
Gene Symbol	Drd5
Uniprot ID	DRD5_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 428-478 on mouse Dopamine Receptor D5 protein.
Immunogen Region	428-478
Specificity	
Immunogen Sequence	