

Dopamine D5 Receptor protein (STJ503946)

STJ503946

GENERAL INFORMATION

Product Type	Proteins
Short Description	Dopamine D5 Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1816
Gene Symbol	DRD5
Uniprot ID	DRD5_HUMAN
Immunogen	Synthetic peptide corresponding to unique sequence on Dopamine Receptor D5 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	