

DDX6 protein (STJ503949)

STJ503949

GENERAL INFORMATION

Product Type	Proteins
Short Description	DDX6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1656
Gene Symbol	DDX6
Uniprot ID	DDX6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 433-483 on probable ATP-dependent RNA helicase DDX6 protein.
Immunogen Region	433-483
Specificity	
Immunogen Sequence	