

DEK protein (STJ503950)

STJ503950

GENERAL INFORMATION

Product Type	Proteins
Short Description	DEK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7913
Gene Symbol	DEK
Uniprot ID	DEK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-330 on C terminal region of Human Draxin protein. Region: NNKCFDDCMC VEGLRCYAKF HRNRRVTRRK
Immunogen	300-330
Region	
Specificity	
Immunogen	
Sequence	