

EDIL3 protein (STJ503951)

STJ503951

GENERAL INFORMATION

Product Type	Proteins
Short Description	EDIL3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13612
Gene Symbol	Edil3
Uniprot ID	EDIL3_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 30-60 on Del-1 protein.
Immunogen Region	30-60
Specificity	
Immunogen Sequence	