

DNA PKcs protein (STJ503959)

STJ503959

GENERAL INFORMATION

Product Type	Proteins
Short Description	DPKcs-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5591
Gene Symbol	PRKDC
Uniprot ID	PRKDC_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 4000-4128 of human DNPK1 protein.
Immunogen Region	4000-4128
Specificity	
Immunogen Sequence	