

C1orf187 protein (STJ503961)

STJ503961

GENERAL INFORMATION

Product Type	Proteins
Short Description	C1orf187-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	374946
Gene Symbol	DRAXIN
Uniprot ID	DRAXL_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-330 on C terminal region of Human Draxin protein. Region: NNKCFDDCMC VEGLRCYAKF HRNRRVTRRK
Immunogen Region	300-330
Specificity	
Immunogen Sequence	