

DRP1 protein (STJ503962)

STJ503962

GENERAL INFORMATION

Product Type	Proteins
Short Description	DRP1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	74006
Gene Symbol	Dnm1l
Uniprot ID	DNM1L_MOUSE
Immunogen	Synthetic peptide within 450-490 on human DRP1 protein
Immunogen Region	450-490
Specificity	
Immunogen Sequence	