

CRMP3 protein (STJ503963)

STJ503963

GENERAL INFORMATION

Product Type	Proteins
Short Description	CRMP3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10570
Gene Symbol	DPYSL4
Uniprot ID	DPYL4_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids 509-527 of Mouse CRMP3. The peptide was covalently modified to achieve the desired antigenicity before coupling to carrier protein. Sequence: PAKPGSGTQARASC SGKIS
Immunogen Region	509-527
Specificity	
Immunogen Sequence	