

Dopamine Receptor D1 protein (STJ503965)

STJ503965

GENERAL INFORMATION

Product Type	Proteins
Short Description	Dopamine Receptor D1-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Drd1
Uniprot ID	DRD1_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	