

## Dopamine Receptor D4 protein (STJ503968)

STJ503968

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Dopamine Receptor D4-protein if for use in research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	WB: 1:1, 000-1:1, 5000
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">13491</a>
<b>Gene Symbol</b>	<a href="#">Drd4</a>
<b>Uniprot ID</b>	<a href="#">DRD4_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	