

Dopamine D5 Receptor protein (STJ503969)

STJ503969

GENERAL INFORMATION

Product Type	Proteins
Short Description	Dopamine D5 Receptor-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:1, 000-1:1, 5000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1816
Gene Symbol	DRD5
Uniprot ID	DRD5_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	