

E2F1 protein (STJ503984)

STJ503984

GENERAL INFORMATION

Product Type	Proteins
Short Description	E2F1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1869
Gene Symbol	E2F1
Uniprot ID	E2F1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids 58-93 of Human E2F1. Sequence: PCDPDLLLFATPQAPRPTPSAPRPALGRPPVKRRLD
Immunogen Region	58-93
Specificity	
Immunogen Sequence	