

E2F3 protein (STJ503988)
STJ503988**GENERAL INFORMATION**

Product Type	Proteins
Short Description	E2F3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13557
Gene Symbol	E2F3
Uniprot ID	E2F3_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 70-120 on E2F3 protein.
Immunogen Region	70-120
Specificity	
Immunogen Sequence	