

E2F4 protein (STJ503989)

STJ503989

GENERAL INFORMATION

Product Type	Proteins
Short Description	E2F4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	104394
Gene Symbol	E2f4
Uniprot ID	E2F4_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 360-410 on mouse E2F4 protein.
Immunogen Region	360-410
Specificity	
Immunogen Sequence	