

E2F5 protein (STJ503990)

STJ503990

GENERAL INFORMATION

Product Type	Proteins
Short Description	E2F5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13559
Gene Symbol	E2f5
Uniprot ID	E2F5_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 5-55 on E2F5 protein.
Immunogen Region	5-55
Specificity	
Immunogen Sequence	