

EHD1 protein (STJ503996)

STJ503996

GENERAL INFORMATION

Product Type	Proteins
Short Description	EHD1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10938
Gene Symbol	EHD1
Uniprot ID	EHD1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 515-565 on human EH-domain containing protein 1.
Immunogen Region	515-565
Specificity	
Immunogen Sequence	