

EHD3 protein (STJ503998)

STJ503998

GENERAL INFORMATION

Product Type	Proteins
Short Description	EHD3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	30845
Gene Symbol	EHD3
Uniprot ID	EHD3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 525-555 on human EH-domain containing protein 3.
Immunogen Region	525-555
Specificity	
Immunogen Sequence	