

EHD4 protein (STJ503999)

STJ503999

GENERAL INFORMATION

Product Type	Proteins
Short Description	EHD4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	30844
Gene Symbol	EHD4
Uniprot ID	EHD4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 420-450 on EH-domain containing 4 protein.
Immunogen Region	420-450
Specificity	
Immunogen Sequence	